

RESTORATIVE



PARADISE
DENTAL TECHNOLOGIES

The PDT Difference

Linda Miller started PDT in Missoula, MT in 2000. It is now the second-largest manufacturer of dental instruments in North America and available in 57 countries.

PDT's resin handles are lightweight and wide with an exaggerated knurled grip, which reduces pinch force and hand fatigue and increases tactile sensitivity.

PDT's instruments are made with durable steel and have stronger, harder tips that stay sharp longer and are easily resharpended.



Proud to Manufacture in Missoula, MT

As you enter our facility in Missoula, MT, you will see a sign that reads: "Beyond these doors, the world's finest dental instrument craftspeople apply their trade."

We love what we do, and we enjoy working with you.

Stay Healthy and Joyful,
Linda Miller
PDT CEO/Founder


rence



hygiene | diagnostic | restorative | surgical | endodontic

ERGONOMIC. EFFECTIVE. EFFICIENT.

DISCOID-CLEOIDS



UWD5
R268




89/92 DC
R271



3/6 DC
R269



4/5 DC R270 Is designed to carve structural features and remove extra material during restorative procedures. This instrument features highly polished tips to prevent material from sticking during procedures.



4/5 DC ☀️
R270



CARVERS



Hollenback 1/2 3 R284 features a round tip designed for contouring occlusal and interproximal composite and amalgam restorations. The highly polished tips allow you to be confident that materials will not stick to the instrument during use.



Tanner 3
R291



IPC Off Angle - Non-Coated
R282



Carver Hollenback 1/2 3 ☀
R284



Ward 1
R292



IPC- Non-Coated
R283



Carver Hollenback 3
R285



Ward 2
R293

BURNISHERS



Burnisher 18
R329



26/27S Burnisher
R332



Acorn 21B ☀️
R335



2/29 Burnisher
R330




27/29 Burnisher
R333



Calcium Hydroxide Placement ☀️
R336



Hollenback 6 ☀️
R331



Ladmore 3
R334



Norton CP
R337



Hollenback 6 R331 is a restorative instrument designed to condense and smooth composite and amalgam material.



Acorn 21B R335 is a burnisher designed with a distinctive acorn shape for the placement and contouring of composite material.



Calcium Hydroxide Placement R336 is designed for the placement of calcium hydroxide or glass ionomer liner in restorative or endodontic dentistry.




PLUGGERS and CORD PACKERS



Blacks Plugger 0/1

R347

Also available in Serrated R347S



Blacks Plugger 1/2

R348

Also available in Serrated R348S



Backaction Plugger 11/12

R349

Also available in Serrated R349S



113 Straight Cord Packer

R363



113 Straight Serrated Cord Packer

R363S



113 Off Angle Cord Packer

R362



113 Off Angle Serrated Cord Packer

R362S

MIXING and EXCAVATORS



38-39 Excavator
R386



63/64 Excavator (15-8-12)
R395



Mixing Spatula/Applicator
R375



#24 Spatula
R374

MARGIN TRIMMERS



28 Margin Trimmer
(10-95-7-14)
R403



29 Margin Trimmer
(10-80-7-14)
R405

CHISELS



11/12 Bi-angle Chisel
(15-8-8)
R419



19 Bi-Bevel Hoe
(3-2-28)
R420

HATCHETS




8-9 Hatchet (10-7-14)
R416



51/52 Hatchett (15-8-12)
R421




COMPOSITE




Gregg 4-5
R294

The gold standard for composite sculpting instruments.


Material will not stick and is easy to manipulate due to titanium nitride coating.



G2 Thin Composite
R286



G3 Thin Composite
R287



G4 Thin Mini Composite
R288

COMPOSITE



CSS-2 Gold Line™ ☀️
R514



CSS-3 Gold Line™ ☀️
R518



CSS-4 Gold Line™
R522



CSS-5 Gold Line™
R526



CSS-6 Gold Line™
R530



COMPOSITE



IPC-T Gold Line™ ☀️
R534



IPC-L Gold Line™ ☀️
R538



8A Gold Line™
R542



CSS-2 Gold Line™ R514 This instrument's hatchet/hoe blades are designed for applying, shaping, and cord packing.



CSS-3 Gold Line™ R518 This instrument's hatchet and plugger are designed for shaping, applying, and condensing.




IPC-L Gold Line™ R538 has a long, thin and flexible blade perfect for composite placement and interproximal contouring.



IPC-T Gold Line™ R534 This thin flexible blade is perfect for composite placement and interproximal contouring.


COMPOSITE



DrV™ 1-2 ☀️
R560



DrV™ 5-6
R562



DrV™ 7-8
R563



DrV™ 3-4 ☀️
R561



DrV™ 1-2 R560 features both tips used for shaping occlusal anatomy in the composite, with one end having a conical shape that can also be used to create the facial anatomy of the upper central incisors.



DrV™ 3-4 R561 features wide and thin hatchet and hoe blades for applying and compacting composite, excellent for broad surfaces.

WAXING



PKT 1 - Thomas Waxing Instrument
R571



PKT 2 - Thomas Waxing Instrument
R572



PKT 3 - Thomas Waxing Instrument
R573



PKT 4 - Thomas Waxing Instrument
R574



PKT 5 - Thomas Waxing Instrument
R575



CROWN REMOVAL



Crown Remover, Vertical
T830




Crown Remover, Horizontal
T831



10g Adhesive Powder
T783



Crown Removing Pliers
T780



Crown Removing Replacement Tips, 2 ea
Blue **Green**
T782 T779



Crown Removing Pliers Kit
T781

RUBBERDAM ACCESSORIES



Rubberdam Punch Forceps
AINSWORTH
T784



Rubberdam Frame
T786




Rubberdam Clamp Forceps
IVORY
T785



RESTORATIVE KITS

Cement Kit R981



Cementing restorations is easy when you have everything in one convenient kit, including a signature Montana Jack® for removing residual cement.

FlipTop™ Cassette 'D' - Fixed Rack 9
T009D-F

Montana Jack® Rigid R138R

Explorer 5 R166

Mirror Handle (Cone-Socket) R245

Mixing Spatula/Applicator R375

Tweezers COLLEGE T040


Articulating Paper Forceps MILLER
T046

#5 PrecisionReflect HD™ Mirror T050

Aspirating Anesthetic Syringe Medium
T235

Bur Block T319

Composite and Crown Prep Kit R984



Designed to complete everyday procedures with ease.

FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider T015E-F

Explorer 5 R166

Mirror Handle (Cone-Socket) R245

IPC- Non Coated R283

G2 Thin Composite R286

G3 Thin Composite R287

113 Straight Serrated Cord Packer R363S

Mixing Spatula/Applicator R375

38-39 Excavator R386

Tweezers COLLEGE T040

Articulating Paper Forceps MILLER T046

#5 PrecisionReflect HD™ Mirror T050


Aspirating Anesthetic Syringe Medium T235

Bur Block T319

IRIS Scissors, curved T800

RESTORATIVE KITS

Crown and Bridge Kit R980



Including PDT's signature Montana Jack® Rigid, the Crown and Bridge Kit streamlines restorative procedures, saving time in the operatory.

FlipTop™ Cassette 'E' - Fixed Rack 15 with
Standard Divider T015E-F

Montana Jack® Rigid R138R

Explorer 5 R166

Probe WRS ACE™-UNC12 R215

Mirror Handle (Cone-Socket) R245

G2 Thin Composite R286

113 Straight Serrated Cord Packer R363S

Mixing Spatula/Applicator R375

38-39 Excavator R386

Tweezers COLLEGE T040

Articulating Paper Forceps MILLER T046


#5 PrecisionReflect HD™ Mirror T050

Aspirating Anesthetic Syringe Medium
T235

Bur Block T319

IRIS Scissors, curved T800

Total Restorative Kit R937



Designed to complete restorative procedures with ease, including placement of amalgam or composite materials.

FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider T015E-F

Montana Jack® R138

Explorer 23/12 Probe R223

Mirror Handle (Cone-Socket) R245

3/6 DC R269

Hollenback 1/2 3 Carver R284

G2 Thin Composite R286

G3 Thin Composite R287

27/29 Burnisher R333

Calcium Hydroxide Placement R336

Blacks Plugger 0/1 R347

113 Off Angle Cord Packer R362

Mixing Spatula/Applicator R375

38-39 Excavator R386

IPC-T Gold Line™ R534

Tweezers MERIAM T042

Articulating Paper Forceps MILLER T046

#5 PrecisionReflect HD™ Mirror T050

Aspirating Anesthetic Syringe Medium T235

Iris Curved Scissors (11.5cm) T800

ANCILLARY RESTORATIVE PRODUCTS



**Amalgam Carrier SS DE
Regular Large**
T351



Amalgam Well
T356



Articulating Paper Forceps MILLER
T046



**Tofflemire Matrix Band
Retainer**
T787

HOME OF THE MONTANA JACK®



THANK YOU FROM ALL OF US AT PDT!



PDT, Inc. is an ISO 13485:2016 registered company.



www.pdtdental.com

(406)626-4153

customer.care@pdtdental.com

8275 US Highway 10 W
Missoula, MT 59808

Ergonomic. Effective. Efficient.