


SURGICAL




PARADISE
DENTAL TECHNOLOGIES

The PDT Differe

Linda Miller started PDT in Missoula, MT in 2000. It is now the second-largest manufacturer of dental instruments in North America and available in 57 countries.

PDT's resin handles are lightweight and wide with an exaggerated knurled grip, which reduces pinch force and hand fatigue while increasing tactile sensitivity.

PDT's instruments are made using durable steel resulting in stronger, harder tips that stay sharp longer, and are resharpended easily.



Proud to Manufacture in Missoula, MT

As you enter our facility in Missoula, MT, you will see a sign that reads: "Beyond these doors, the world's finest dental instrument craftspeople apply their trade."

We love what we do, and we enjoy working with you.

Stay Healthy and Joyful,
Linda Miller
PDT CEO/Founder

erence



hygiene | diagnostic | restorative | surgical | endodontic

ERGONOMIC. EFFECTIVE. EFFICIENT.

PAGE 2


SCALPEL HANDLES



Scalpel Handle

T704

A high-quality stainless steel scalpel handle made for use with standard scalpel blades.



Scalpel Handle, Adjustable in 7 Positions

T705

Adjusts to 7 different angles for better access to hard-to-reach areas. Easily disassembles for cleaning and sterilization.

ELEVATORS

APEXO



APEXO Elevator 302
T115



APEXO Elevator 303
T116



APEXO Elevator 301
T117

CRYER



Root Elevator CRYER 39 Smooth Handle ☀️
T108S



Root Elevator CRYER 40 Smooth Handle
T109S

ELEVATORS

SELDIN



**Root Elevator Seldin No. 301 Straight
Smooth Handle**
T102S



Root Elevator 46R serrated
T105



Root Elevator Seldin 34, Smooth Handle
T103S



Root Elevator Seldin 34S, Smooth Handle
T104S

LUXATING INSTRUMENTS



Luxating Instrument - 2mm
T120



Luxating Instrument - 3mm ☀️
T121



Luxating Instrument - 4mm
T122



Luxating Instrument - 5mm
T123

PDT's Luxating Instruments feature a fine cutting edge that gently severs the periodontal ligament during tooth extraction. Not to be used as an elevator.

ROOT TIP PICKS

HEIDBRINK



Root Elevator HEIDBRINK L Curved ☀️
T110



Root Elevator HEIDBRINK R Curved ☀️
T111




Root Elevator HEIDBRINK Straight
T112

HEIDBRINK




Elevator Fragment Ejector/Root Tip Pick
HEIDBRINK
T320

PERIOSTEAL ELEVATORS



2mm Periosteal Elevator
R442



3mm Periosteal Elevator
R443



Periosteal Elevator PRICHARD
T740



Periosteal Elevator Molt 9, hollow handle
T118

PERIOTOMES

Micro-Serrated Edges to Assist in Extractions

*Not to be used as an elevator.



Micro-Serrated Periotome™ - Anterior
R256



Micro-Serrated Periotome™ - Posterior
R258



Why PDT Periotomes?

- Less damage to alveolar bone and better clinical results.
- Vertically insert the micro-serrated edge to cut, not tear, the periodontal fibers.
- Minimize post-surgical bone loss and tissue trauma.
- Less traumatic for the patient and faster healing.

EXTRACTING FORCEPS



Extracting Forceps 150 Extra Grip ☀️
T790



Extracting Forceps 150S Extra Grip ☀️
T791



Extracting Forceps 151 Extra Grip ☀️
T792



Extracting Forceps 151S Extra Grip ☀️
T793



Extracting Forceps 23 Cowhorn ☀️
T794



Extracting Forceps 69
T798

EXTRACTING & ROOT TIP FORCEPS



Extracting Forceps 88L ☀️
T795



Root Forceps TC ☀️
T797



Extracting Forceps 88R ☀️
T796



Stieglitz 30° Splinter Forcep
T732

RONGEURS & MALLETS



Bone Rongeur Friedman
T730



Mallet - 180g
T869



Bone Rongeur Micro-Friedman ☀️
T731

SURGICAL CURETTES



The Slade Blade™ was designed in collaboration with periodontist Scott Froum, DDS.


SLADE BLADE™




Slade Blade™ 1 Anterior Narrow
R432



Slade Blade™ 2 Anterior Wide
R433



Slade Blade™ 3 Posterior Narrow
R434



Slade Blade™ 4 Posterior Wide
R435



Slade Blade™ 5 Posterior Wide
R436

LUCAS



Lucas 86
R437

NEEDLE HOLDERS



**TC Needle Holder Castroviejo, Straight 17.5
T088**



**TC Needle Holder Castroviejo Straight 14
T710**



**TC Needle Holder CRILE WOOD
T714**



**TC Needle Holder MATHIEU
T715**

SCISSORS

METZENBAUM




Scissors METZENBAUM Pointed, Curved
T808

GOLDMAN-FOX



Scissors Goldman-Fox Curved 12.5
T813



Scissors Ligature 13
T816

BEEBEE



BEEBEE Crown Scissors, Curved
T804

SCISSORS

IRIS



IRIS Scissors, Curved ☀️
T800

LA-GRANGE



Scissors La-Grange 11.5
T820



IRIS Scissors, Curved Black
T801



Scissors La-Grange 11.5 cm Super Cut ☀️
T821



IRIS Scissors, Straight
T802

HEMOSTATS & TISSUE FORCEPS

HEMOSTAT



Halsted-Mosquito Curved 12
T716

ADSON



Tissue Pliers ADSON
T080



Tissue Pliers ADSON (with tooth)
T081



Spoon Pliers Hemingway
T862

TISSUE FORCEPS

SEMKEN-TAYLOR



Tissue Pliers SEMKEN-TAYLOR 15
T082



Tissue Pliers SEMKEN-TAYLOR 15 (with tooth)
T083

CORN SUTURE



Sutures Pliers Corn
T084



Micro Surgical Corn Suture Plier ☀️
Diamond Dusted
T085

MICRO STRAIGHT



Pliers Micro Straight 17.5 cm
T086

MICRO STRAIGHT



Pliers Micro Straight 17.5 (Diamond Coated)
T087

PERIODONTAL CHISELS & KNIVES




#13/13KL Surgical Chisel
R446



End Cutting Intrasulcular Knife
R444




Chisel OCHSENBEIN 1
T775



Modified Orban Knife
R445



Chisel OCHSENBEIN 2
T776



19/20 USC Towner Knife
R449

KITS

Basic Extraction Kit

R983



Lucas 86

R437

Tweezers COLLEGE

T040

#5 PrecisionReflect™ HD Mirror

T050

Mirror Handle (Cone-Socket)

R245

Tissue Pliers ADSON

T080

Retractor Caewood Minnesota 16

T100

Root Elevator Seldin No. 301 Straight

T102S

Root Elevator Seldin 34S

T104S

Periosteal Elevator Molt 9 Hollow Handle

T118

Aspirating Anesthetic Syringe Medium

T235

**Elevator Fragment Ejector/Root Tip Pick,
HEIDBRINK**

T320

Scalpel Handle

T704

TC Needle Holder CRILE WOOD

T714

Scissors Goldman-Fox Curved 12.5

T813

Miller-Colburn Bone File #1, Serrated

T850

FlipTop™ Cassette 'E' - Utility

T015E-U


Silicone Finger Mat

T018FM

KITS

Tissue Grafting Kit

R927



FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider T015E-F

Mirror Handle (Cone-Socket) R245

#5 PrecisionReflect™ HD Mirror T050

Younger-Good 7-8 R077

Explorer 23/Probe UNC-15 R226

Suture Pliers Corn T084

Pliers Micro Straight 17.5 T086

TC Needle Holder Castroviejo, Straight 17.5 T088

Adjustable Scalpel Handle T705

Scissors Goldman-Fox curved 12.5 T813

2mm Periosteal Elevator R442

3mm Periosteal Elevator R443

End Cutting Intrasulcular Knife R444

Modified Orban Knife R445

TISSUE PUNCHES



4mm

5mm

6mm

PDT Tissue Punches are available in three sizes: 4mm, 5mm, and 6mm. Each size is designed to precisely harvest tissue grafts, tailored to the specific area of the mouth. Using the appropriate size ensures optimal fit and coverage at the recipient site, while minimizing unnecessary tissue removal from the donor area for efficient and effective procedures.



Tissue Punch, Keyes 4mm
T859



Tissue Punch, Keyes 5mm
T860



Tissue Punch, Keyes 6mm
T861

ASPIRATING ANESTHETIC SYRINGES



Aspirating Anesthetic Syringe (Small)
T234



Aspirating Anesthetic Syringe (Medium)
T235



Aspirating Anesthetic Syringe (Large)
T236

Stainless Steel and Made in the USA



small

medium

large

The PDT Aspirating Anesthetic Syringe is designed to allow you to aspirate comfortably with an ergonomic thumb ring and finger rest (in three sizes).

Improved carpule visibility due to viewing windows on both sides of barrel.

Lightweight for better control and comfort.

BONE CARRIERS



Bone Carrier plugger 2mm w/ depth marking
T856



Bone Carrier plugger 3mm w/ safety stop
T857



Bone Carrier plugger 4mm w/ safety stop
T858



Carrier BONE - 3.5 mm
T845



Carrier BONE - 4.5 mm
T846



Carrier BONE - 6.0 mm
T847

BONE GRAFTING ACCESSORIES



Mixing Bowl for bone material
T849



Measuring/mixing bowl 25cc
T848



Bone Scraper
T844



Miller-Colburn File #1, Serrated
T850



Bone Mill with Titanium beaks
T843

RETRACTORS



Retractor Caewood Minnesota 16
T100



Retractor VESTIBULUM 15
T101



Proudly woman-owned, family-run and American made!



PDT, Inc. is an ISO 13485:2016 registered company.

www.pdtdental.com

(406)626-4153

customer.care@pdtdental.com

8275 US Highway 10 W
Missoula, MT 59808

Ergonomic. Effective. Efficient.