

# KITS



**PARADISE**  
DENTAL TECHNOLOGIES

# The PDT Difference

Linda Miller started PDT in Missoula, MT in 2000. It is now the second-largest manufacturer of dental instruments in North America and available in 57 countries.

PDT's resin handles are lightweight and wide with exaggerated knurling, which reduces pinch force and hand fatigue and increases tactile sensitivity.

PDT's instruments are made with durable steel and have stronger, harder tips that stay sharp longer and are easily resharpened.




## Proud to Manufacture in Missoula, MT

As you enter our facility in Missoula, MT, you will see a sign that reads: "Beyond these doors, the world's finest dental instrument craftspeople apply their trade."

We love what we do, and we enjoy working with you.

Stay Healthy and Joyful,  
Linda Miller  
PDT CEO/Founder

# rence



hygiene | diagnostic | restorative | surgical | endodontic

**ERGONOMIC. EFFECTIVE. EFFICIENT.**

# RDH EXAM KIT

US CDT Code D0120/0150

**RDH Exam Kit**  
R972



**FlipTop™ Cassette 'A' - Fixed Rack 5**  
T005A-F

**Probe WRS ACE™ - UNC 12** R215

**Explorer 5** R166

**Explorer ODU 11-12** R167

**Mirror Handle (Cone-Socket)** R245

**#5 PrecisionReflect™ HD Mirror** T055

A separate exam kit allows you to:

- Select the correct scalers to meet the needs of each patient
- Save money - fewer diagnostic instruments needed

# PEDIATRIC PROPHY

US CDT Code D1120

## Pediatric Prophylaxis R973



**FlipTop™ Cassette 'A' - Fixed Rack 5**  
T005A-F

**204S** R135

**Columbia 13-14** R064

**Jack B. Nimble®** R130

**N128-L5 Mini** R093

Pair with the RDH Exam Kit for a full suite of diagnostic and hand instruments. The Pediatric Prophylaxis Kit has everything a clinician needs for primary and mixed dentition for smaller pediatric or early adolescent dentition.

# ADULT PROPHY

US CDT Code D1110

**Adult Prophylaxis**  
R974



**FlipTop™ Cassette 'A' - Fixed Rack 5**  
T005A-F

**Montana Jack®** R138

**204S** R135

**Amazing Gracey™ 15-16** R030

**Amazing Gracey™ 13-14** R028

**Barnhart 5-6** R061

Featuring three signature instruments and two clinician favorites, this kit makes adult or adolescent hygiene procedures a breeze.

# SCALING & ROOT PLANING

US CDT Code D4341/4342

## SCALING & ROOT PLANING R975



FlipTop™ Cassette 'B' - Fixed Rack 7 T007B-F

Queen of Hearts™ R144

Montana Jack® Pink R138BC

Montana 4-5™ R120

Amazing Gracey™ 15-16 Extended Reach R044

Amazing Gracey™ 13-14 Extended Reach R043

Columbia 4L-4R R067

N128-L5 Mini R093

Tweezers COLLEGE T040

When treating periodontal disease, instruments with unique designs and longer shanks are invaluable. With a 6mm-long cutting edge and ideal angulation, the Queen of Hearts™ is designed for hard-to-reach subgingival surfaces and furcations. The Montana 4-5™, with its extended, rigid shank, provides greater interproximal adaptability for the hard-to-reach molars.

# PERIO MAINTENANCE

US CDT Code D4910

## Perio Maintenance R976



**FlipTop™ Cassette 'B' - Fixed Rack 7**  
T007B-F

**Queen of Hearts™** R144

**Montana Jack® Pink** R138BC

**23 Slim (Pointed Toe)** R098S

**Amazing Gracey™ 15-16**  
**Extended Reach Mini** R008

**Amazing Gracey™ 13-14**  
**Extended Reach Mini** R007

**Barnhart 5-6** R061

**N128-L5 Mini** R093

It is essential to have the right instruments for Periodontal Maintenance. This kit features clinician favorites with extended shanks and mini cutting edges for hard-to-reach subgingival surfaces.



# GINGIVAL INFLAMMATION

US CDT Code D4346

## Gingival Inflammation R977



**FlipTop™ Cassette 'A' - Fixed Rack 5**  
T005A-F

**Montana Jack® Pink** R138BC

**N2** R139

**Columbia 13-14** R064

**N1** R096

**Jack B. Nimble® Slim** R130S

Featuring both signature instruments and clinician favorites, like the Montana Jack® and the Columbia 13-14, this kit has everything a clinician needs to treat both healthy or inflamed gingiva.

This kit is also great for orthodontics and lingual bars.

# GROSS DEBRIDEMENT

US CDT Code D4355

## Gross Debridement R978



**FlipTop™ Cassette 'B' - Fixed Rack 7**  
T007B-F

**Montana Jack® Rigid** R138R

**Montana 4-5™** R120

**Amazing Gracey™ 15-16 Rigid** R030R

**Amazing Gracey™ 13-14 Rigid** R028R

**Columbia 4L-4R** R067

Effortlessly remove tenacious calculus utilizing instruments featuring extended, rigid shanks and signature designs, like the Montana 4-5™ and Amazing Graceys™.

# I DREAM OF HYGIENE™

## I DREAM OF HYGIENE™ KIT R913



Are you looking for a kit that does it all? The I Dream of Hygiene™ Kit showcases top-ranked diagnostic, sickle and curette designs to treat any patient.

**FlipTop™ Cassette 'C' -  
Double Rack 18 T009C-D**

**Barnhart 5-6 R061**

**Langer 3-4 R080**

**Jack B. Nimble® Slim R130S**

**Montana Jack® R138**

**Queen of Hearts™ R144**

**Explorer ODU 11-12 R167**

**Explorer 23/Probe 12 R223**

**Mirror Handle (Cone-Socket) R245**

**#5 PrecisionReflect™ HD Mirror T050**

**One-Piece College Plier with U Bend  
T043**

# UNIVERSAL SRP KIT C

## Universal SRP Kit C R915



FlipTop™ Cassette 'B' - Double Rack 14 T007B-D

Columbia 4L-4R R067

Langer 3-4 R080

H5-L5 R090

Montana Jack® R138

Explorer ODU 11-12 R167

Probe ACE™- 12 Short R216

Mirror Handle (Cone-Socket) R245

#5 PrecisionReflect™ HD Mirror T050

Tweezers COLLEGE T040

Smaller kits can pack a punch!

The Universal SRP C Kit features top clinician-ranked curettes and sickles in addition to diagnostic instruments to treat a healthy adult or an adult with gingival inflammation.

Double rack cassette allows clinicians to add more of their favorites.

# JUST BE NICE KIT

## JBN KIT R970



**FlipTop™ Cassette 'B' - Fixed Rack 7 T007B-F**

**Montana Jack® R138**

**Jack B. Nimble® Slim R130S**

**Amazing Gracey™ 11-12 Extended Reach Mini R006**

**Amazing Gracey™ 13-14 Extended Reach Mini R007**

**Columbia 13-14 R064**

**Explorer 23/Probe OW R224**

**Mirror Handle (Cone-Socket) R245**

**#5 PrecisionsReflect™ HD Mirror T050**

The Just Be Nice Kit allows for agility throughout any hygiene procedure. Featuring signature sickles and Extended Reach Mini Amazing Graceys™ in addition to a clinician top choice, Columbia 13-14, plus light-weight diagnostic instruments. Also available without the cassette, item number R971.

# PINEYRO ARCH™ KIT

## Pineyro Arch™ Kit R930



Created in  
collaboration with  
Dr. Alfonso Pineyro



CELLERANT  
BEST OF CLASS  
HYGIENE AWARD

2021



**FlipTop™ Cassette 'A' - Fixed Rack 5  
T005A-F**

**Pineyro Arch™ Ti 1 R630**

**Pineyro Arch™ Ti 2 R631**

**Pineyro Arch™ Ti 3 R632**

**Pineyro Arch™ Ti 4 R633**

**Pineyro Arch™ Ti Kit  
(Instruments Only)  
R931**

Removal of biofilm, calculus, debris or cement from implant collars and threads is no longer a challenge with the specially designed Pineyro Arch™ Kit. Save chair time during preventive procedures by eliminating the need to remove a fixed prosthesis at each visit.

# WINGROVE™ TITANIUM IMPLANT INSTRUMENT “GO-TO” SET

## Wingrove™ “Go-To” Set R900



Created in  
collaboration with  
Susan Wingrove, BS,  
RDH



FlipTop™ Cassette ‘A’ - Fixed Rack 5  
T005A-F

Wingrove™ B5-6 Ti R661

Wingrove™ L3-4 Ti R680

Wingrove™ N128 Blunt-Langer 5 Mini Ti  
R693

DiamondHead™ Ceramic Stone T061

EasyView™ Probe UNC12 T216

Efficiently assess and maintain any  
implant with the Wingrove™ “Go-To”  
Set.

# DDS EXAM KIT

**DDS Exam Kit**  
R979



**FlipTop™ Cassette 'A' - Fixed Rack 5**  
T005A-F

**Explorer 5** R166

**Probe WRS ACE™-UNC 12** R215

**Tweezers COLLEGE** T040

**#5 PrecisionReflect™ HD Mirror** T050


**Mirror Handle (Cone-Socket)** R245

Complete comprehensive, periodic  
or limited exams with ease.



# CEMENT KIT

## Cement Kit R981



Cementing restorations is easy when you have everything in one convenient kit, including a signature Montana Jack® for removing residual cement.

**FlipTop™ Cassette 'D' - Fixed Rack 9**  
T009D-F

**Montana Jack® Rigid** R138R

**Explorer 5** R166

**Mirror Handle (Cone-Socket)** R245

**Mixing Spatula/Applicator** R375

**Tweezers COLLEGE** T040

**Articulating Paper Forceps MILLER**  
T046


**#5 PrecisionReflect™ HD Mirror** T050

**Aspirating Anesthetic Syringe Medium**  
T235

**Bur Block** T319

# COMPOSITE AND CROWN PREP

## Composite and Crown Prep Kit R984



Designed to complete everyday procedures with ease.

**FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider** T015E-F

**Explorer 5** R166

**Mirror Handle (Cone-Socket)** R245

**IPC- Non Coated** R283

**G2 Thin Composite** R286

**G3 Thin Composite** R287

**113 Straight Serrated Cord Packer** R363S

**Mixing Spatula/Applicator** R375

**38-39 Excavator** R386

**Tweezers COLLEGE** T040

**Articulating Paper Forceps MILLER** T046

**#5 PrecisionReflect™ HD Mirror** T050


**Aspirating Anesthetic Syringe Medium** T235

**Bur Block** T319

**IRIS Scissors, curved** T800

# CROWN AND BRIDGE KIT

## Crown and Bridge Kit R980



Including PDT's signature Montana Jack® Rigid, the Crown and Bridge Kit streamlines restorative procedures, saving time in the operatory.

FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider T015E-F

Montana Jack® Rigid R138R

Explorer 5 R166

Probe WRS ACE™-UNC12 R215

Mirror Handle (Cone-Socket) R245

G2 Thin Composite R286

113 Straight Serrated Cord Packer R363S

Mixing Spatula/Applicator R375

38-39 Excavator R386

Tweezers COLLEGE T040

Articulating Paper Forceps MILLER T046

#5 PrecisionReflect™ HD Mirror T050


Aspirating Anesthetic Syringe Medium T235

Bur Block T319

IRIS Scissors, curved T800

# TOTAL RESTORATIVE KIT

## Total Restorative Kit R937



Designed to complete restorative procedures with ease, including placement of amalgam or composite materials.

FlipTop™ Cassette 'E' - Fixed Rack 15  
with Standard Divider T015E-F

Montana Jack® R138

Explorer 23/Probe 12 R223

Mirror Handle (Cone-Socket) R245

3/6 DC R269

Carver Hollenback 1/2 3 R284

G2 Thin Composite R286

G3 Thin Composite R287

27/29 Burnisher R333

Blacks Plugger 0/1 R347

113 Off Angle Cord Packer R362

Mixing Spatula/Applicator R375

38-39 Excavator R386

IPC-T Gold Line™ R534

Tweezers MERIAM T042

Articulating Paper Forceps MILLER T046

#5 PrecisionReflect™ HD Mirror T050

Aspirating Anesthetic Syringe Medium T235

Calcium Hydroxide Placement R336

Iris Curved Scissors (11.5cm) T800

# BASIC EXTRACTION KIT

## Basic Extraction Kit R983




Streamline extraction procedures by having the most frequently utilized instruments at your fingertips.

FlipTop™ Cassette 'E' - Utility T015E-U  
Tissue Pliers ADSON T080  
Retractor Caewood Minnesota 16 T100  
Root Elevator Seldin No. 301 Straight T102S  
Root Elevator Seldin 34S T104S  
Periosteal Elevator Molt 9 Hollow Handle T118  
Lucas 86 R437  
Tweezers COLLEGE T040  
Silicone Finger Mat T018FM  
Mirror Handle (Cone-Socket) R245  
#5 PrecisionReflect™ HD Mirror T050  
Aspirating Anesthetic Syringe Medium T235  
Elevator Fragment Ejector/Root Tip Pick, HEIDBRINK T320  
Scalpel Handle T704  
TC Needle Holder CRILE WOOD T714  
Scissors Goldman-Fox Curved 12.5 T813  
Miller-Colburn Bone File #1, Serrated T850

# ENDO KIT

## Endo Kit R982



Streamline endodontic procedures by having the most frequently utilized instruments at your fingertips.

**FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider T015E-F**

**Mirror Handle (Cone-Socket) R245**

**DG 16 Explorer R311**

**D11 Endo Spreader R312**

**1/3 Endo Plugger R314**

**5/7 Endo Plugger R315**

**9/11 Endo Plugger R316**

**Glick 1 R317**

**31L Endo Excavator R318**

**Locking Cotton Plier T045**

**Articulating paper forceps MILLER T046**


**#5 PrecisionReflect™ HD Mirror T050**

**Aspirating Anesthetic Syringe Medium T235**

**Rubberdam clamp forceps IVORY T785**

# TISSUE GRAFTING KIT

## Tissue Grafting Kit R927



Streamline tissue grafting procedures by having the most frequently utilized instruments at your fingertips.

**FlipTop™ Cassette 'E' - Fixed Rack 15 with Standard Divider** T015E-F

**Mirror Handle (Cone-Socket)** R245

**#5 PrecisionReflect™ HD Mirror** T050

**Younger-Good 7-8** R077

**Explorer 23/Probe UNC-15** R226

**Suture Pliers Corn** T084

**Pliers Micro Straight 17.5** T086

**TC Needle Holder Castroviejo, Straight 17.5** T088

**Adjustable Scalpel Handle** T705

**Scissors Goldman-Fox curved 12.5** T813

**2mm Periosteal Elevator** R442

**3mm Periosteal Elevator** R443

**End Cutting Intrasulcular Knife** R444

**Modified Orban Knife** R445

# GLEASON GUIDE™ SHARPENING SYSTEM



The Gleason Guide™ was created in collaboration with Roger Gleason.



**Gleason Guide™**  
T065

Stainless steel instrument sharpening guide.


**Make sharpening easy: stop, rest, glide!**

Three Beveled edge guides: one side Gracey, one Universal/Sickle and the other to shape the toe and sharpen O'Hehirs, Chisels and Hoes.



**Diamond Head Stone**  
T061

Fine ceramic stone that will not pit or groove. No need to use oil or water to sharpen. This stone is designed to maintain the shape and edge of your instruments.



**Transformation Stone**  
T062

Coarse ceramic stone that will not pit or groove. No need to use oil or water. This stone is designed to reshape the edge of your instruments.



# ULTIMATE EDGE™ KITS



**Ultimate Edge™  
Sharpening Kit**  
T066



**Ultimate Edge™ Kit  
with Transformation Stone**  
T068




**Ultimate Edge™  
Transformation Kit**  
T067

Save time and increase accuracy with the Gleason Guide™. Sharpen any instrument using the three beveled edges.

All components, including the PingRing™, are autoclavable.

# PDT'S INSTRUMENTS ARE *Anatomically Color Coded* for faster recognition and setups



# HOME OF THE MONTANA JACK®



*Proudly woman-owned, family-run and American made!*

PDT, Inc. is an ISO 13485:2016 registered company.





[www.pdtdental.com](http://www.pdtdental.com)

(406)626-4153

[customer.care@pdtdental.com](mailto:customer.care@pdtdental.com)

8275 US Highway 10 W  
Missoula, MT 59808

**Ergonomic. Effective. Efficient.**